



New DNA High Performance Air Filter Release Sheet • #25 / 2013

Dear Business Partners and Friends,

We proudly announce the official release of the *P-SY2SC13-01* 

DNA air filter for the SYM models:

- ✓ JOYRIDE 125 03'-08'
- / JOYRIDE 200 03'-07'
- / JOYRIDE 200 EFI 07'-10'
- ✓ JOYRIDE 200 EFI EVO 09'-13'
- This filter features DNA<sup>®</sup>'s advanced FCd<sup>1</sup> (Full Contour design).
- A laser cut supporting frame is included with the DNA filter.

• Installation of this new DNA Air filter is very easy; simply follow the installation instructions included in the SYM workshop manual.

• The filtering efficiency<sup>2</sup> is extremely high at **98-99%** filtering efficiency (ISO 5011), with 4 layers of DNA<sup>®</sup> Cotton.

• The flow of this new DNA Fcd filter is high, +59.41% more

than the SYM stock paper filter!

DNA Fcd air filter flow: **149.50 CFM** (Cubic feet per minute) @1,5"H2O corrected @ 25degrees Celsius SYM stock paper filter: **93.78 CFM** (Cubic feet per minute) @1,5"H2O corrected @ 25degrees Celsius


• This DNA<sup>®</sup> filter is designed as a High flow Air filter for: **'Road & Off road use'**.





FCd (Full Contour design) is the innovative design by DNA<sup>®</sup>, that allows the filtering material to follow precisely the contour of the air box and uses the complete air box

surface as "active filtering area" eliminating "dead spots" that cause turbulence, increasing air flow and filtering efficiency.















AIR FLOW DATA MEASURED WITH DNA'S ROTRONICS FLOWSCAN COMPUTERIZED FLOWBENCH

