

Dear Business Partners and Friends,

We proudly announce the official release of the **P-SY1SC22-01** DNA High Performance Air Filter for the following

SYM models:

- JET X 125 21'-22'
- JET 14 125/200 17'-22'
- SYMPHONY 50/125/200 15'-22'



- This filter features DNA®'s advanced FCd<sup>1</sup> (Full Contour design).
- Installation of this DNA Air filter is very easy; simply follow the installation instructions included in the SYM workshop manual.
- A perfect airtight sealing and trouble free filter installation for the user is guaranteed by using a high quality **EVA seal**.
- The filtering efficiency<sup>2</sup> is extremely high at **98-99%** filtering efficiency (ISO 5011), with 4 layers of DNA® Cotton.
- DNA Filters are supplied pre-oiled from our factory. The correct amount of DNA oil is applied on each filter individually.
- The flow of this new DNA filter is **+22.49%** more than the stock filter!  
- DNA air filter's flow: **76.80 CFM**  
(Cubic feet per minute) @1,5"H<sub>2</sub>O corrected @ 25degrees Celsius
- - SYM stock filter's flow: **62.70 CFM**  
(Cubic feet per minute) @1,5"H<sub>2</sub>O corrected @ 25degrees Celsius
- This DNA® filters is designed as High Flow Air filter for: **'Road and City use'**.


**SYM**

**JET X / JET 14 / SYMPHONY MODELS**

DNA PART N°:

**P-SY1SC22-01**



---

STOCK FILTER'S AIR FLOW

**62.70 CFM** ✓

---

DNA® FILTER'S AIR FLOW

**76.80 CFM** ✓

---

DNA® INCREASED AIR FLOW

**+22.49 %** ✓


---


DNA® FILTERING EFFICIENCY

**98-99%** ✓

---

EAN No : **5212008312046**

**1**  FCd (Full Contour design) is the innovative design by DNA®, that allows the filtering material to follow precisely the contour of the air box and uses the complete air box surface as "active filtering area" eliminating "dead spots" that cause turbulence, increasing air flow and filtering efficiency.

**2**  Filtering efficiency is the amount of "dirt" the filter can maintain (stop) and protect the engine efficiently. For example the DNA® Filter for every 100 grams of dirt that it will receive, it will hold **98-99 grams**, this applies even to fine dirt as small as 5 microns.

